


Welcome to the latest improvements in VetConnect* PLUS

Save the tests that matter most to you

You can save tests as you order them in the IDEXX Diagnostics window or in the Directory of Services.



To find and retrieve a saved test, click **Saved**.

Click the star ★ to save a test. You can save as many tests as you'd like.

Review other recent improvements below

Print or email results from the Home page or the search results page

Print or email a single result

The screenshot shows a 'Complete' status page with a 'MOST RECENTLY UPDATED' header. It lists three test results:

- BENTLEY COLE**: 19979999 Senior Profile with Triglycerides-St...
9999 SAMPLE/TEST INFO NEEDED
Sep 14, 2017
- MOLLY ALLEN**: 16019999 Total Health™ PI
Sep 14, 2017
- CAESAR REMINGTON**: 16019999 Total Health™ Plus with Triglyceri...
Sep 14, 2017

On the right side, there are buttons for 'Generate Diagnostic Summary', 'IDEXX Reference Laboratory Operational Alerts', and 'See what's new on VetConnect PLUS'. A red arrow points to the action icons (print, email, mark as read) for Molly Allen's result.

Hover the cursor over a test result to display print, email, and mark-as-read buttons.

Print or email multiple results

The screenshot shows a table of search results with a toolbar at the bottom. The table has columns for patient name, ID, veterinarian, request number, and date. The toolbar includes buttons for 'Print', 'Email', 'Mark Read', and 'Deselect All'. A red arrow points to the 'Mark Read' button.

Patient Name	ID	Veterinarian	Request	Date	Action
KANEL HONTORIA	582-1	OSCAR PAEZ CASTRO	Req: 118	Apr 10, 2013 12:29 pm	<input checked="" type="checkbox"/>
DUTCH PEREZ	598-1	OSCAR PAEZ CASTRO	Req: 141	May 2, 2013 7:00 pm	<input checked="" type="checkbox"/>
MATHEW DIAZ	634-1	Shih Tzu	Req: 198	Jul 13, 2013 10:40 am	<input type="checkbox"/>
TOKY MONTES	84-1	Shih Tzu	Req: 221	Jul 30, 2013 4:20 pm	<input type="checkbox"/>
POLO 18339 TERRY 15128	Yorkshire Terrier	RUBY		Nov 23, 2012 6:55 pm	<input type="checkbox"/>
1	Maltese			Sep 3, 2012 1:14 pm	<input type="checkbox"/>

Select up to ten results, and then print, email, or mark as read with one click.

Filter results by diagnostic testing type

The screenshot shows the IDEXX VetConnect PLUS interface. At the top, there are navigation links for 'Home', 'Directory of Services', 'Imaging', and 'Telemedicine'. Below that, there's a 'DIAGNOSTICS' section with a search bar and a date range '9/7/2017 - 9/14/2017'. There are three filter buttons: 'Reference Laboratory', 'In-house Laboratory', and 'Imaging'. A red circle highlights the 'Imaging' button, and a red arrow points to it from the text below. Below the filters, there are two columns: 'In Process' and 'Complete', each with a 'MOST RECENTLY UPDATED' header and a list of results.

Select an option to filter results by test type.

Track test progress at a glance

In Process

MOST RECENTLY UPDATED

LILLY COOLIDGE 5580-1
4033 Urine Culture and MIC Susceptibility*
Nov 26, 2015
Ordered

MOLLY BOTT 16038-1
VetLyte Electrolyte Analyzer
ProCyte Dx Hematology Analyzer
Catalyst Dx Chemistry Analyzer
Nov 26, 2015
Running

DIVA CHESTER 34910-3
401 Aerobic & Anaerobic Cultures*
Nov 26, 2015
At the Lab

CHARM FOSTER 21277-1
4033 Urine Culture and MIC Susceptibility*
Nov 26, 2015
Partial

- ✓ Hematology
- ✓ Chemistry
- ✓ Urinalysis (updated)
- ☒ Endocrinology (updated)

New status icons keep you up to date on test status.

Icon...	Applies to...	Indicates....
Ordered	IDEXX Reference Laboratories	Order has been created
At the Lab	IDEXX Reference Laboratories	Order received but no results available yet
Partial	IDEXX Reference Laboratories	Some tests complete; partial results available
Running	IDEXX in-house analyzers	Test is running now; no results available yet

To view partial test results:

Place the mouse pointer over **Partial** to see which test categories are complete and which are pending. Click any category, partial or complete, to see available test results.

CHARM FOSTER 21277-1
4033 Urine Culture and MIC Susceptibility*
Nov 26, 2015
Partial

- ✓ Hematology
- ✓ Chemistry
- ✓ Urinalysis (updated)
- ☒ Endocrinology (updated)


In this example, three categories are complete and Endocrinology is pending.

Easily find and add tests to samples already at IDEXX Reference Laboratories

Initial results often indicate more testing is needed. We've made it quick and easy to search all IDEXX tests and add those you need.


1. Open the patient's test results in VetConnect PLUS.
2. Click the new **Add to Order** button at the top. The button is active only if add-on tests are available for this order.
3. Search all IDEXX tests, and then select the ones to add. You'll receive a confirmation message.
4. (Optional) Print the updated test requisition form for your records.

You **don't** need to send the printed form to IDEXX. Your test request is sent electronically, and we already have your samples.



The screenshot shows the IDEXX VetConnect PLUS interface for a patient named MACY GILKS. At the top, there is a navigation bar with the date 2017 APR 22 and a blue button labeled '+ Add to Order'. Below this, there is a table of test results. The table has three columns: Test Name, Value, and Reference Range. The tests listed are RBC, Hematocrit, and Hemoglobin. The RBC value is 7.6, Hematocrit is 0.46, and Hemoglobin is 177. The reference ranges are 5.4 - 8.7 x10¹²/L for RBC, 0.38 - 0.57 L/L for Hematocrit, and 134 - 207 g/L for Hemoglobin. The interface also shows a date 4/22/17 (Order Received) and a timestamp Yesterday 6:16 AM (Last Updated).

Test Name	Value	Reference Range
RBC	7.6	5.4 - 8.7 x10 ¹² /L
Hematocrit	0.46	0.38 - 0.57 L/L
Hemoglobin	177	134 - 207 g/L



The screenshot shows a confirmation message in the IDEXX VetConnect PLUS interface. The message is displayed in a white box with a yellow checkmark icon. The text reads: 'You have requested additional testing to reference laboratory samples'. Below this, it states: 'Your request will be processed by IDEXX Customer Support. If you have any questions or need assistance please call, 1-800-667-3411.' At the bottom of the message, it says: 'If you have already submitted samples, you do not need to send a new copy of your requisition form.' At the bottom of the interface, there are two buttons: 'Done' and 'Print Requisition'.